# ADACIRCUIT

# ADACIRCUIT

## 1. Intended use and general warnings

- This product is only used for extracorporeal blood circulation during Adacolumn<sup>(1)</sup> granulocyte and monocyte apheresis with the Adamonitor SC. Any other use must be avoided.
- This product has been sterilized by Ethylene Oxide
- This device is a disposable product and should not be re-used or re-sterilized. The reuse of the device may pose risks related to microbiological and/or physico-chemical aspects.

## 2. Precautions before use

- The Adacircuit package must be checked for damage before use; if the package is damaged in any way the Adacircuit should not be used.
- In order to avoid contamination the Adacircuit should be used immediately after removal from the packaging.

## 3. Safety during set-up

### Cautions

- When setting up the blood circuits with the blood pump care should be taken to avoid any contamination.
- The blood tubings are distinguished by colours. The blood inflow line has red caps and connectors and the blood outflow line has blue caps and connectors.
- Make sure to connect each line of the Adacircuit correctly to the Adacolumn, the Adamonitor SC and the patient.
- For safety reasons the Adacircuit should be primed by pumping approximately 1000 mL saline through the tubings before use.

### Procedure

*Note: Please refer to the Adacolumn package insert and Adamonitor SC manual details*

- Attach the pump tube section of the Adacircuit (thick section of tubing which is part of the blood inflow line) to the blood pump.
- Prime the Adacircuit by pumping approximately 1000 mL of saline through the tubings.
- When priming with saline is completed, heparinised saline should be pumped through the Adacircuit to replace normal saline. The drip chamber of the unit should be filled up to approximately 80% (four fifths).
- Check that there is no leakage in any part of the Adacircuit.
- If anything abnormal is noticed with the Adacircuit during the priming procedure, the circuit should be discarded and not used.

Please note: Execution of the procedure is the responsibility of the qualified medical professional performing the procedure. Description of the instruction for use is given as information only. The physician must evaluate the appropriate procedure based on their knowledge and experience. The instructions for use of the blood pump and the Adacolumn must be read before the procedure is started.

## 4. Care when starting the procedure

- Adacircuit should be used only under the supervision of a qualified medical professional.
- When using the pump tube, do not conduct circulation above 300 mL/min for more than 8 hours.
- When clamping the tube be careful not to damage the circuit.
- If blood leaks out and/or air bubbles are seen in the tubings, the circulation should be discontinued and the Adacircuit exchanged.

- If coagulation occurs, circulation should be discontinued and corrective actions should be initiated at the discretion of the physician (e.g. exchanging the circuit and restarting the apheresis).

## 5. Care after use

- After the apheresis has been completed, sufficient re-infusion of residual blood should be performed.
- Once re-infusion is completed, the circuit should be removed from the patient, disconnected from the equipment, and disposed appropriately as medical waste.
- When disposing the used product, be careful to avoid contamination.

## 6. Precautions for storage and handling

- Do not use if the packing or the device is wet.
- The device is to be stored in a clean dry place.
- Avoid exposure to direct sunlight, ultraviolet rays, high temperatures or humidity.

## 7. Warranty Limits

The Adacircuit has been carefully designed, manufactured, and tested according to the manufacturer's specifications;

However, it is possible that for unexpected reasons the product may not perform as designed.

The information in the attached package insert is detailed and is an integral part of the declaration of warranty.

If the product is used in any other way other than instructed by Bioengineering Laboratories in this package insert, Bioengineering Laboratories does not accept any liability.

The above paragraph is not intended to contradict local laws where the product is distributed.

In the event that any part of the above paragraph is judged to be in conflict with local laws, the validity of the other sections are not compromised and all rights and obligations will be taken and applied as if the above instructions do not include the conflicting part.

### Manufacturer:

BIOENGINEERING LABORATORIES SRL

Via A. Vivaldi 32/A,

22063 Cantù (CO) - ITALY

Phone: +39 0317377737

Fax: +39 031730467 [www.bioengineeringlab.com](http://www.bioengineeringlab.com)

<sup>(1)</sup> Adacolumn is a trade mark of JIMRO Co. Ltd. Japan

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